

Clariant Produkte (Deutschland) GmbH Analytical Sciences, G860 Dipl.-Ing. (FH) Martina Hirschen

Tel. 069 305 16861 Date : 2025-07-16 Page : 1/2

REPORT OF ANALYSIS - PHYSICS, PARTICLE & SURFACE ANALYSIS

Sample: TB1180B

Order title/Product: Graphene / Raman microscopy

Project: Graphene

Account: Curia Germany GmbH
Customer: Dr. Thorsten Busch
Order no.: 25-050404-2

Remarks: Raman measurement on request of the customer

Method:

Raman spectroscopy/microscopy using Senterra I Raman microscope (Bruker).

Procedure:

Principle:

Objective: Raman spectroscopy

Field of Application: Powder

The substance is applied to a quarz microscope slide and examined under the

microscope. The objective is selected and focused to examine an exact location

using Raman. At this position, several points are examined with the laser and a

corresponding number of Raman spectra are recorded.

Reference: N/A

Experimental Details and Result

Equipment: Senterra I Raman microscope (Bruker)

Magnification: 20 x

Wavenumber range: $45 - 3700 \text{ cm}^{-1}$

Wavenumber step size 3 – 5 cm⁻¹

Co-addition: 100

Exposure time: 10000 ms

Evaluation: Averaging over similar spectra

Electronic document valid without signature!



Clariant Produkte (Deutschland) GmbH Analytical Sciences, G860 Dipl.-Ing. (FH) Martina Hirschen Tel. 069 305 16861

Date : 2025-07-16 Page : 2/2

Result:

Absorption band	2D band	G band	D band	D/G ratio
Integration range	2791 - 2557	1753 – 1474	1474 – 959	
Integral Typ A	9,690	226,6	322,3	1,422
Peak high Typ J	0,123	1,895	1,708	0,901

Spectra see attachment.




Figure 1: TB1180B, average spectrum over 8 recordings